



New Patents Pending

Dental Treatment Face Mask



Inventor
Y.C. Cheng
(Mask / マスク)

A whole new dental treatment experience, right in your face!



Structure

- Inner layer :
absorbent micro-fibered
non-woven
- Middle layer :
super fluid absorbent layer
- Outer layer :
protective non-woven

Types :

- (1) Regular
- (2) Advanced
- (3) Activated Carbon
- (4) Eye Protection Bundled

Categories :

- (1) Adult - (a) Regular (M) (b) Extended (L)
- (2) Child (S)

Features :

- (1) Cutting-edge Innovation
Mouthpiece opening facilitates dental treatment and scaling, while maintaining nasal protection.
- (2) Face fit & protection
Evenly distributed mount provides complete facial cover against instruments and pollutants.
- (3) Stain-resistant & infection-free
Super fluid absorbent layer draws out polluted fluid, while effectively prevents contamination.